

Product details

Information for disposable stretchers sheet cover

- 1) Color: white,blue,green,etc
- 2) Size: 30"x72",30"x84",31"x77", 34"x90",can be made according to client's demand
- 3) Style: fitted corner
- 4) Free samples of non woven bed sheet available for test, samples can be offered in 3~5 days.

Fabric of cot sheet

1. PP

Fabric weight : 20-65gsm.

PP material is made from hydrophobic polypropylene .

Latex-free;abrasion-resistant; low lint;with a high level of fluid repellency.

2. SMS

Fabric weight : 20-65gsm.

SMS is made from hydrophobic SMS/Spun lace material.

Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to alcohol ,blood, body fluids and pathogens.

3) Microporous

Fabric weight : 45-65gsm.

The material is waterproof and breathable,good permeability, with high resistance to impact mechanical properties of the toughness, texture.

4) PP+PE

Material weight:40-65gsm.

Polyethylene film is coated in polypropylene non-woven fabric.

Latex-free; abrasion-resistant; completely impervious fluid and alcohol repellent.

The trade terms of disposable bed sheets

MOQ	5,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	10,000pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	According to the quantity, normally 30days after deposit

Product display



The regular shipment way



Our business process

