

Chinese best suppliers for medical surgical earloop face mask , nonwoven face mask , breathing filter mask

Specifications of medical surgical earloop face mask

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	White/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing design	All inner box and outer carton printing can be as your design
Applications	Hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on
MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

Packaging of nonwoven face mask

50pcs/box,40boxes/ctn,52x38x30cm;

1*20 GP container can hold 1,000,000 Pcs of masks;

1*40 HC container can hold 2,450,000 Pcs of mask.

We can also do packaging of breathing filter mask as your requirements.

Discarded

- 1) Once finished use can be recycled or incinerated.
- 2) Packaging should be deposited in the appropriate container for recycling

Picture of breathing filter mask



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

