

Introduction of disposable clip cap:

Brand Name: ORIENT HONEST

Model Number: OH-1001, OH1002, OH1003

Certification: CE / ISO13485 / FDA

Place of Origin: China

Material : PP nonwoven

Color: white,green,blue,pink,yellow,pink,orange,purple,red (all colors can be available)

Style: elastic clip caps

Weight: 10g-30g

Size: 18''□19''□21''□etc

Packing details: 100 pcs/bag,20 pcs/ctn

Specifications of disposable clip cap:

Product:	medical consumables high quality colorful non-woven single or double elastic 19" 20" 21" custom size disposable clip cap
Material:	Nonwoven pp with elastic
Grammage:	10g~16g
Style:	Single elastic or double elastic
Size:	19", or 18",20",21", etc.
Color:	White, Green, yellow, Blue, etc.
MOQ	100,000pcs
Packaging:	100pcs/bag ; 2000pcs/ctn
Application:	for household, medical, surgical, pharmaceutical, laboratorial, industrial sectors, painting, cleaning room, etc.
Storage:	Stored in dry, humidity below 80%, ventilated, non-corrosive gases warehouse
Terms of payment:	T/T, L/C, Paypal,D/P,D/A,etc.
Other service:	OEM and Customized Logo/brand
Certification:	CE, ISO13485,FDA

Features of disposable clip cap:

1. Breathable, non-absorbent, durability
2. With single or double elastic available, in pleated or clip shape with various colours available.
3. Bacteria filtration efficiency (BFE) has reached 99.1%. This rate is much higher than normal standard of 97%.

4. Free samples and customization available
5. Clip-shaped cap before unfolded
6. We could guarantee best quality, fast delivery, development new material.

Pictures of disposable clip cap:





orient Honest
MEDICAL INTERNATIONAL





