

The property of Cotton facial towel

All cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:


1. Size: 20 x 20 cm each piece or customized
2. Made of cotton material, soft and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .
4. Can be wet and dry, are great for fragile or sensitive skin
5. Great water absorption, skin-friendly and comfortable.
6. Cotton tissue can instead of face towel and will be much cleaner to dry your face after washing.

The Specification of Cotton facial towel

Name	Cotton facial towel
Material	100% cotton; spunlace nonwoven
Style	Plain, mesh, pearl and customization patterns available
Size	15*20cm, 20*20cm ,or other size can be customized
Color	white
Packing	1 roll/portable bag ;
Carton Box	All inner bag and carton printing can be as your design
Carton Size	52x38x30cm /52x19x30cm
Applications	personal care, makeup removing, house cleaning, travel use, camping, similar outdoor activities and business occasions.
MOQ	1000 bags
Sample	Can be supplied free for your quality checking within 3 working days
Payment way	T/T, PAYPAL, Westunion, L/C, D/A, D/P, ESCROW
Delivery time	as you quantity
FOB Port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of Cotton facial towel





COTTON TISSUE

- No chemical addition
- No fiber falls on skin
- No any harm to the skin



orient Honest
MEDICAL INTERNATIONAL



COTTON TISSUE

No chemical addition
No fiber falls on skin
No any harm to the skin




Plant fibres



Safe enough


7.8inch



Elastic



Skin friendly



Fibrous porous structure



The power of touch



 Remove the bacteria



 Clean computer screen and accessories



 Clean cellphone screen



 Clean Automobile Interiors



orient Honest
MEDICAL INTERNATIONAL

Double-side design:

Side A:

Stronger & More Durable



Side B:

Silky feeling for better
care of your skin

