

The property of Disposable Black medical mask


All masks are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Super soft earloop makes comfortable wearing for long time

The Specification of dDisposable Black medical mask

<u>Fabric Material</u>	Layer 1: PP non woven Layer 2: high filtration paper Layer3: pp non woven
Style	3 ply ,with round ear loop
Dimensions size	17.5*9.5 cm
Fabric weight	25gsm + 25gsm + 25gsm, could be as your request
Color	white, baby blue, baby pink, grey
Certificate	CE, ISO13485, FDA
MOQ:	200000pcs
Packing design	can be as your design
Packing	50 pcs/box,2000 pcs/ctn, or as your requirement
Capability	2 million pcs per day
Delivery time	7-30 days
<u>Applications</u>	hospital, laboratory, clinic, food industry, dust-free workshop

The picture of Disposable Black medical mask



Disposable black surgical mask



- 1. 3 layers : PP + filter paper + PP
- 2. BFE >99%
- 3. EN 14683, CE FDA ISO 13485
- 4. Size : 17.5*9.5cm







