

Nonwoven Disposable Caps, Disposabel Nurse Cap, Surgical Bouffant Cap With Single Elastic

It is made of PP non woven , sewing and automatic style\

Advantages:

- 1, Good protective, light weight, latex free.
- 2,Sewing style, Comfortable wearing.
- 3,White, Blue, Green, Yellow, Pink, Orange, Red colors can be available.

Indications

For hospital clinical medical and patients, dust free work shop, food factory.

Clip caps packing

100PCS/bag, 2000pcs/carton. Carton size: 58*25*27cm 1*20ft can load 715CTNS, total 1430000 pcs. 1*40HQ can load 1737CTNS, total 3474000pcs.

Name	Disposable surgical non woven elastic caps
Material	PP 10-30g
Style	Bouffant cap, clip cap, mob cap, doctor cap
Color	white, green, blue, yellow, pink, red, orange, etc.
Packing	100pcs/bag, 2000pcs/ctn
packing design	all inner bag and carton printing can be as your design
MOQ	100,000pcs
Sample	can be supplied free for your quality checking within 3 days.
Payment way	T/T, Paypal, Western union, L/C,D/A, D/P
Delivery time	Capacity :200,000pcs/day, 1*40HQ can be finished within 10-15days.
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of disposable bouffant cap















<http://https://membercenter.made-in-china.com/>