

#### Key Specifications/Special Features:

- 1.Product name: disposable facial towel
- 2.Shelf life: 3 years
- 3.Size information: 100x200mm,150x200mm,200x200mm,other size can be customized
- 4.Specification: 100 pcs/ bag
- 5.Main Ingredient: 100% cotton
- 6.Application: Makeup remover, Baby Skin care, Daily care...

#### Key features:

- 1.100% cotton,can add hot water to keep baby warm in cold winter
- 2.Without any chemical substances
- 3.Also suitable for adults to clean face and hands

#### Packing:

- 1.Cardboard box, can add inner transparent polybag if needed
- 2.Sheet size: 100x200mm,150x200mm,200x200mm

#### Our services:

- 1.Professional designer team offer you high-standard & satisfactory OEM packages design proposal
- 2.OEM services: weight, color, material, styles can be customized
- 3.Arrange shipment for you with competitive freight charge
- 4.Brand marketing promote your sale

The complete closed type, dust-proof and bacteria-free production workshop with central air condition guarantees all the products made in the international standard environment

#### Specification

Product Name	disposable facial towel
Material	100% cotton
Application	Makeup remover, Baby Skin care, Daily care
Feature	Natural,Organic,Lint-free,Safe, Eco-friendly, 100% biodegradable and 100% compostable
MOQ	1000Packs
Packaging	Customization
BasisWeight	30gsm-120gsm
Sheet Size	(100mm-200mm)*200mm
Pattern	Plain, mesh, pearl and customization patterns available