

## Key Specifications/Special Features:

### 1. PRODUCT SIZE:

S M L XL XXL

### 2. MATERIAL:

PP 30-70GSM

SMS 30-70GSM

SF 50-70CM

PP+PE 30-60GSM

### 3. MATERIAL TREATMENT:

Material TREATMENTRESULT

ANTI-WATER:YES

ANTI-DUST:YES

ANTI-STATIC:YES/NO

ANTI-BLOOD:YES/NO

ANTI-ALCOHOL:YES/NO

### 4. COLOR: WHITE,BLUE,DARK BLUE,ORANGE AND SO ON

### 5. USEAGE:

A. KIMONObe used in the spa and sauna to convenience the body against the water and the dust.

B. KIMONO can be used in the massage to relax your body.

C.KIMONO can be used in the cleaning room,for example the GMP.

D.KIMONO can be used in the hospital.

### 6. PACKING:

10pc/bag ,100pcs/ctn

### 7. OEM SERVICE:

A.Can print the band for our customer on the PE bag

B.Can print the band for our customer on the Carton box

### 8. SAMPLES:

Samples can be sent out with in 3days after you confirmed.

### 9. DELIVERY DATE:

It can be sent in 15days after received the deposite.

product name:	Disposable Kimono
Size:	S:110*130cm;M:115*137cm,L:120*140cm,XL:125*145cm;XXL:130*150cm Or customized
Color:	white,blue,pink,yellow or customized

Material:	SBPP,SMS,Spunlace
Sample:	Free samples available
Application:	Mainly used in SPA ,beauty salon hospital etc
MOQ:	5000PCS
Package:	10pc/bag ,100pcs/ctn
Payment:	T/T,L/C
Delivery time:	Shipping in 20 days