

All lab coat are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

Fabric Material	SMS
Style	With ties on collar and waist, knitted cuffs
Dimensions size	S-3XL
Fabric weight	35-65g, could be as your request
Color	Blue color
Certificate	ISO13485, FDA
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	10-30 days
Applications	laboratory, clinic, food industry, dust-free workshop,hospital

The picture of Disposable Lab coat with knitted cuffs medical manufacturer




orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL







orient Honest
MEDICAL INTERNATIONAL

