

All medical gowns are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

Fabric Material	SMS material
Style	With ties on collar and waist, knitted cuffs
Dimensions size	120*140 cm ,acn be customized
Fabric weight	35-45g, could be as your request
Color	Blue color
Certificate	ISO13485, FDA,CE
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	7-30 days based on the order quantity
Applications	laboratory, clinic, food industry, dust-free workshop,hospital

The picture of Disposable Level 1 medical gown with knitted cuffs medical supplier



Tie on neck



Tie on belly



Tie on waist




Knitted cuff





**Velcro
back neck**



Round collar



Tie on waist



Knitted cuff

