

## **Nonwoven Disposable Caps, Disposable Nurse Cap, Surgical Bouffant Cap With Single Elastic**

It is made of PP non woven, sewing and automatic style.

### **Advantages**

1. Good protective, light weight, latex free.
2. Sewing style, Comfortable wearing.
3. White, blue, green, yellow, pink, orange, red colors can be available.

### **Indications**

For hospital clinical medical and patients, dust free work shop, food factory.

### **Clip caps packing**

100PCS/bag, 2000PCS/carton. Carton size: 58X25X27cm 1\*20 container can load 715ctns, total 1430, 000PCS. 1\*40 container can load 1737ctns, total 3474, 000PCS

Name Disposable surgical non woven elastic mob caps

Material PP 10-30g

Style Bouffant cap, clip cap, mob cap. mop cap. doctor cap

Size 18",19",20",21",22",24"

Color white,green,blue,yellow,pink,red,orange etc.

packing 100 pcs/bag,2000pcs/ctn packing design

all inner bag and carton printing can be as your design

MOQ 100,000 pcs Sample can be supplied free for your quality checking within 3 days Payment way

T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,

Delivery time Capacity: 200,000pcs/day, 1\*40HQ can be finished within 10-15 days;

FOB port Wuhan/Shanghai/Ningbo/Guangzhou