

Nonwoven Disposable Caps, Disposable Nurse Cap, Surgical Bouffant Cap With Single Elastic

It is made of PP non woven, sewing and automatic style.

Advantages

1. Good protective, light weight, latex free.
2. Sewing style, Comfortable wearing.
3. White, blue, green, yellow, pink, orange, red colors can be available.

Indications

For hospital clinical medical and patients, dust free work shop, food factory.

Clip caps packing

100PCS/bag, 2000PCS/carton. Carton size: 58X25X27cm 1*20 container can load 715ctns, total 1430, 000PCS. 1*40 container can load 1737ctns, total 3474, 000PCS

Name Disposable surgical non woven elastic mob caps

Material PP 10-30g

Style Bouffant cap, clip cap, mob cap. mop cap. doctor cap

Size 18",19",20",21",22",24"

Color white,green,blue,yellow,pink,red,orange etc.

packing 100 pcs/bag,2000pcs/ctn packing design

all inner bag and carton printing can be as your design

MOQ 100,000 pcs Sample can be supplied free for your quality checking within 3 days Payment way T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,

Delivery time Capacity: 200,000pcs/day, 1*40HQ can be finished within 10-15 days;

FOB port Wuhan/Shanghai/Ningbo/Guangzhou