

Indications of disposable shoe covers:

For hospital clinical medical and patients ,dust free work shop ,food factory.

Packing of disposable shoe covers:

100pcs/bag,2000pcs/ ctn. carton size:36\*29\*54cm.,1\*40container can load 240,0000pcs.

Name	shoe cover
Material	PP,PP+PE,CPE
Style	Sewing ,automatic machine made
Size	38x15cm,40x15cm,41x17cm.....
Color	White,blue,green,red colors can be available.
packing	100pcs/bag,2000pcs/ctn
packing design	all inner box and carton printing can be as your design
MOQ	50000 pcs
Sample	can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 600,000pcs/day, 1*40HQ can be finished within 30 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou





orient Honest  
MEDICAL INTERNATIONAL





orient Honest  
MEDICAL INTERNATIONAL



