

Disposable Protective [Non Woven Lab Coat Black Color Suppliers](#)

- 1, Free sample.<?XML:NAMESPACE PREFIX = O />
- 2, Small order acceptable;
- 3, L/C, T/T, D/P terms available;
4. FOB, CNF, CIF, EXW terms available;
- 5, Competitive and stable price because we are factory.
- 6, Qualified products. 3rd part product inspection is acceptable.
- 7, OEM service. More options on material, pattern, size and packing;
- 8, Timely delivery. The best lead time is 25-30 day for a 40HQ container.
- 9, Professional export service. Documentation, clearance, sea transportation, all in one.
- 10, Good after sale service. We take care complaints as we would like to be your long-term supplier!

Specification of [Non Woven Lab Coat Black Color](#)

Name	Disposable Protective Non Woven Lab Coat Black Color Suppliers
Size	S-5XL sizes are available, regular size. S: 98x118cm, M: 103x123cm,L: 107x128cm, XL: 110x134cm, XXL: 117x140cm. (Making production according to client's samples & specifications)
Material	PP45g
Style	Stand collar, Zip front, Elastic cuffs
Color	black
Application	Medical & health / Household / Laboratory
Characteristic	Durable, Eco-friendly, Non-toxic, Breathable and Flexible
	Anti-alcohol, Anti-Blood, Anti-Oil
	Waterproof, Acid proof, Alkali proof
MOQ	5000 PCS
Sample	can provide within 2days
Payment term	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Delivery time	One 40hq container can be finished within 20-30days.



orient Honest
MEDICAL INTERNATIONAL

