

### 1. PP Protection non sterile surgical gowns

It is made from hydrophobic polypropylene material, Latex-free; Abrasion-resistant; Low lint; With a high level of fluid repellency.

Color: White, green, blue, yellow, orange and so on

Material weight: 30-65GSM.

### 2. PP+PE protection non sterile surgical gowns

It is made from **SMS** material, Latex-free; Abrasion-resistant; Completely impervious fluid and alcohol repellent.

Color: White, blue.

Material weight: 40-65GSM.

### 3. SMS protection non sterile surgical gowns (Standard and Reinforced):

It is made from hydrophobic SMS/Spunlace material, Latex-free; Abrasion-resistant; Low lint; With a high level of fluid repellency; A good barrier to blood, body fluids and pathogens.

With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.

Available color: White, blue.

Material weight: 35-65GSM.

Tie on/Velcro on the neck, elastic/cuffed cuff, standard or reinforced.

## The Specification of gown

Name	<b>Disposable SMS surgical gown</b>
Material	SMS, SMMS,
Style	Tie on/Velcro on the neck, elastic/cuffed cuff, standard or reinforced
Size	S-5XL
Color	white,green,blue,yellow,pink,red, grey,black etc.
packing	10pcs/box,10 bags/ctn or 1 pcs/bag,50 pcs/bag
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Sample	can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,
Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou