

Product details

Description of disposable bed sheets for hospital

1. Fabric type: PP,SMS,Microporous,spunlace
2. Fabric weight: 20-100 GSM
3. Color: White, blue, yellow, green or etc...
4. Size: 210X90X20cm, 190X90X15cm or etc...
5. Style: elastic in two ends,or elastic in four corner,we also have flat style
6. Anti-Static treated: available
7. Flame-resistant treated: available, we add the fire-resistant material in the fabric, making the fabric has high resistance to the fire,and the fire can be extinguished in 3-5 seconds.

Description of disposable surgical pillow case

1. Fabric type: PP,SMS,Microporous,Spunlace
2. Fabric weight: 20-100 GSM
3. Color: White, blue, yellow, green or etc...
4. Size: 32X42cm+10cm; 40X70cm+10cm; 60X60cm, or etc...
5. Seal: ultrasonic or sewing seam
6. Style: open style and cover in the opening

The trade terms of hotel bed sheet set

MOQ	5,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	10,000pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	According to the quantity,normally,30 days after deposit

Product display



The regular shipment way



FAQ

Q:1. Which are the main countries you export to?

A: Germany, Switzerland, Austria, Italy, France, UK,USA, Saudi Arabia, Chile, Mexico, Peru, Brazil, Korea, Japan ,Thailand and so on.

Q:2. Could you provide us an overview of the range of the products that you specifically manufacture from your own factory?

A: Our own factory mainly produce all kinds of disposable face masks,caps, coveralls, surgical gowns, lab coats, bed sheets, arm sleeves, shoe covers, surgical drapes, rain ponchos ,aprons and so on.

Q: 3. How many production lines are there for your major range of the products?

A: 2 production lines for fabric PP ,SMS, PP/PE material.

2 dust free workshops for surgical sterile products

5 automatic caps machines and 5 automatic shoe cove machines.

10 automatic face mask machines.