

Specification of Disposable double layers stretcher cover with absorb& water proof:

1.-Material: PP+Microporous

2.-Size: Width 90cm(+5cm), Length 210cm(+5cm)

3.-Color: Blue White Green

4.-Packing: 10pcs/polybag, 100pcs/ctn.

5.-Certification: CE ISO FDA

Function	Absorb& Water proof
Type of Shipping	Sea & air & express shipping
Packing	Plastic bag / paper carton or customized
Carton Size	48x39x10cm or customized
Product Size	70x170cm,80x180cm,90x210cm or customized
Delivery Time	About 7-15 days, based on order quantity
MOQ	10000pcs
Material	<u>PP+Microporous</u>
Sample	Free Sample
Payment Terms	TT/LC/Western union/PayPal
Color	White/Blue/Green
Port of shipment	Shanghai/Ningbo/Wuhan/Shenzhen

Function:

1.-double layered structure with absorb & water proof function

2.-Stretchable & Fixable band should attached

3.-Bursting Strength: Over 125kPa

4.-Tensile strength: (MD 170N), (CD 110N)



orient Honest
MEDICAL INTERNATIONAL



