

## Product Description

Disposable 3 Layers Non Woven Anti Virus Ear Loop Protection High Filtration Face Mask Suppliers made of 3-layers,inside and outside layers are PP non woven,middle layer is high standard 99 filter paper.With Ear-loop style.

## Advantages

- 1.Directly manufacturer,ISO CE FDA certificated 3ply non woven disposable surgical face mask , any colors.
- 2.ISO,CE,FDA, BFE>99%,type II /Type IIR standard.
- 3.Capacity: 1000,000pcs/day, 1\*40HQ can be finished within 10-15 days;
- 4.Laxet free ,easy breath, eco friendly, OEM packing design, service in 24 hours.
- 5.Skin-friendly, Hypoallergenic ,Comfortable, High breathability, fiberglass free, latex free

## Specifications

Material	3ply, PP non woven 25g+filter 25g+pp non wovn 25g
Feature	Skin-friendly, Hypoallergenic ,Comfortable, High breathability, fiberglass free, latex free
Style	Elastic ear-loop (latex free)
Size	17.5x9.5cm(adult)
Color	white,green,blue,pink,yellow,pink,orange,purple,red,apple green,black,(all colors can be available)
packing	Non-sterile:50 pcs/box,2000pcs/ctn or 1000 pcs/ctn
packing design	All inner box and carton printing can be as your design
carton size	52x38x30cm
Applications	hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	10000 PCS
Sample	Can be supplied free for your quality checking within one week
Payment term	T/T,L/C,D/A,D/P/ESCROW
Delivery time	1*40HQ can be finished within 7-10 days.as you quantity( 2 millon pcs per day)
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

## Packaging & Shipping

Non sterile: 50 pcs/box, 40 boxes/ctn or 20 boxes/ctn

## Our Services

### 1. Provide the highest quality products

Professional QC team strictly control each step of the production until loading.

### 2. The fastest delivery time

If urgent, 1\*40HQ can be finished within 5-7days!

Normally, 1\*40HQ can be finished within 10-15 days!

Busy season, 1\*40HQ can be finished within 25-30 days!

### 3. The most reasonable price

Continuously invest in improving our skilled workers, products and processes to minimize the cost of each step!

### 4.The sweet and quick service

Our workers will be always so patiently making a response of your any questions within 5-10 minutes.

### 5.We can also help to source connected products.

Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!

## Company Information

Wuhan Orient Honest Medical Products Co.,ltd ( " Orient Honest") is a global product development and the directly manufacturing company that specializes in innovation and cost effective medical and protective supplies. Orient Honest has expanded to our extensive product lines such as all kinds of disposable protective coveralls, surgical gowns, lab coats, face masks, bed sheets, surgical caps, shoe covers, arm sleeves, rain ponchos, aprons, gauze swabs, cotton rolls and so on.

Orient Honest is a fast-growing company, customer satisfaction is always our mission and goals. Our commitment to continuous improvement of our manufacturing processes, enhancing product performance and developing new cost-effective products to meet the demands of various customers and markets. All of our products comply with international quality standards, such as ISO13485,CE, FDA.

Our goods are already exported to Europe, North America, South America, Asia, and so on.