

## Our Advantage

1. Directly manufacturer
2. ISO13485, CE,FDA certificate
3. Capacity: 10, 000PCS/day, 1\*40HQ can be finished within 10-15 days;
4. OEM acceptable, service in 24 hours.

Disposable surgical gowns are made from high-performance isolation protective material. Widely used in medical safety protective, aseptic workshop, protective isolation, mining electronics factory, food factory farm animal husbandry biohazard and so on.

With two kinds processes: Sewing or ultrasound thermal bonding.

## Material Features

### 1. PP Disposable surgical gowns

It is made from hydrophobic polypropylene material, Latex-free; Abrasion-resistant; Low lint; With a high level of fluid repellency.

Color: White, green, blue, yellow, orange and so on

Material weight: 30-65GSM.

### 2. PP+PE Disposable surgical gowns

It is made from PP+PE material, Latex-free; Abrasion-resistant; Completely impervious fluid and alcohol repellent.

Color: White, green, blue, yellow.

Material weight: 40-65GSM.

### 3. SMS Disposable surgical gowns (Standard and Reinforced):

It is made from hydrophobic SMS/Spunlace material, Latex-free; Abrasion-resistant; Low lint; With a high level of fluid repellency; A good barrier to blood, body fluids and pathogens.

With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.

Available color: White, blue.

Material weight: 35-65GSM.

Tie on/Velcro on the neck, elastic/cutted cuff, standard or reinforced.

Name	Disposable isolation gown
Material	PP, PP/PE, SMS, SMMS, CPE
Style	Tie on/Velcro on the neck, elastic/cutted cuff, standard or reinforced
Size	S-5XL
Color	white,green,blue,yellow,pink,red, grey,black etc.
packing	10pcs/box,10 bags/ctn or 1 pcs/bag,50 pcs/bag
packing design	all inner box and carton printing can be as your design
MOQ	5000 pcs
Sample	can be supplied free for your quality checking within 3 days
Payment way	T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,

Delivery time	Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days;
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou

<http://https://www.china-disposable.com/products/Disposable-surgical-gowns.html#.YljnArUzaUk>