

Disposable polypropylene carrying sheet for repositioning and transferring patients

Features of disposable carrying sheet

- 1) The disposable patient carrying sheet is specially designed from a non-woven, water repellent polypropylene material.
- 2) It's incredibly strong with a maximum user weight of 400kg.
- 3) Each sheet is made to make turning, repositioning and transfer of patients very easy with 6 handles on each side.

The specification for polypropylene carrying sheet

Material	polypropylene
Color	blue
Style	Ultrasonic heat sealing ,with hole in the length side
Size	195cm*85cm (can be made according to client's demand)
Application	Hospital, emergency center, clinic
MOQ	5,000 Pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	10,000pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	According to the quantity, normally 30days after deposit
Packing details	Sterile: 1 Pcs/paper plastic bag,50 Pcs/Ctn Non-sterilize: 1 Pcs/bag, 50 Pcs/ctn Or as your requirements

Pictures of disposable polypropylene carrying sheet



orient Honest
MEDICAL INTERNATIONAL

