

All Surgical gowns are made of the highest quality materials, providing excellent protection and superior comfort:

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

<u>Fabric Material</u>	SMS material
Style	With ties on collar and waist, knitted cuffs
Dimensions size	X-3XL
Fabric weight	30-60g, could be as your request
Color	Blue color
Certificate	ISO13485, FDA
MOQ:	10000 pcs
Packing design	can be as your design
Packing	1pc/plastic paper bag,50pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	7-30 days
Applications	laboratory, clinic, food industry, dust-free workshop,hospital

[The picture of Disposable surgical gown with knitted cuffs medical manufacturer](#)





Tie on neck



Tie on belly



Tie on waist



Knitted cuff

