

**All Surgical gowns are made of the highest quality materials, providing excellent protection and superior comfort:**

- Non-linting materials are smooth and comfortable to eliminate skin irritation
- Highly breathable construction layers remain cool, even during long cases
- Various Style can be customized according to your different requirement

<b>Fabric Material</b>	SMS material
<b>Style</b>	With ties on collar and waist, knitted cuffs
<b>Dimensions size</b>	X-3XL
<b>Fabric weight</b>	30-60g, could be as your request
<b>Color</b>	Blue color
<b>Certificate</b>	ISO13485, FDA
<b>MOQ:</b>	10000 pcs
<b>Packing design</b>	can be as your design
<b>Packing</b>	1pc/plastic paper bag,50pcs/ctn, or as your requirement
<b>Capability</b>	60,000pcs per day
<b>Delivery time</b>	7-30 days
<b>Applications</b>	<a href="#">laboratory, clinic, food industry, dust-free workshop,hospital</a>

**[The picture of Disposable surgical gown with knitted cuffs medical manufacturer](#)**





**Tie on neck**



**Tie on belly**



**Tie on waist**



**Knitted cuff**