

## Our Advantage

1. Directly manufacturer
2. ISO13485, CE,FDA certificate
3. Capacity: 10, 000PCS/day, 1\*40HQ can be finished within 10-15 days;
4. OEM acceptable, service in 24 hours.

Disposable surgical gowns are made from high-performance isolation protective material. Widely used in medical safety protective, aseptic workshop, protective isolation, mining electronics factory, food factory farm animal husbandry biohazard and so on.

With two kinds processes: Sewing or ultrasound thermal bonding.

## Material Features

### 1. PP Disposable surgical gowns

It is made from hydrophobic polypropylene material, Latex-free; Abrasion-resistant; Low lint; With a high level of fluid repellency.

Color: White, green, blue, yellow, orange and so on

Material weight: 30-65GSM.

### 2. PP+PE Disposable surgical gowns

It is made from PP+PE material, Latex-free; Abrasion-resistant; Completely impervious fluid and alcohol repellent.

Color: White, green, blue, yellow.

Material weight: 40-65GSM.

### 3. SMS Disposable surgical gowns (Standard and Reinforced):

It is made from hydrophobic SMS/Spunlace material, Latex-free; Abrasion-resistant; Low lint; With a high level of fluid repellency; A good barrier to blood, body fluids and pathogens.

With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.

Available color: White, blue.

Material weight: 35-65GSM.

Tie on/Velcro on the neck, elastic/cutted cuff, standard or reinforced.

|                |  |
|----------------|--|
| Name           | Disposable surgical gowns  |
| Material       | PP, PP/PE, SMS, SMMS, CPE  |
| Style          | Tie on/Velcro on the neck, elastic/cutted cuff, standard or reinforced |
| Size           | S-5XL  |
| Color          | white,green,blue,yellow,pink,red, grey,black etc.                      |
| packing        | 10pcs/box,10 bags/ctn or 1 pcs/bag,50 pcs/bag                          |
| packing design | all inner box and carton printing can be as your design                |
| MOQ            | 5000 pcs   |
| Sample         | can be supplied free for your quality checking within 3 days           |

|               |  |
|---------------|--|
| Payment way   | T/T, PayPal, Western union, L/C,D/A,D/P,ESCROW,                    |
| Delivery time | Capacity: 10,000pcs/day, 1*40HQ can be finished within 10-15 days; |
| FOB port      | Wuhan/Shanghai/Ningbo/Guangzhou                                    |