

All shoe cover are made of the highest quality materials, providing excellent protection and superior comfort:

- With effective protection.
- Soft non-irritating elastic band.
- fabric for comfort and economy.
- Customers logos can be printed on packings

Faberic Material:	PE, CPE
Style:	Handmade, Made by machine are available
Dimensions size:	15*40 cm, 17*40 cm, 17*42 cm etc.
Color:	Blue , Pink , Green
Certificate:	CE, ISO13485, FDA
MOQ:	100000 PCS
Packing design:	can be as your design
Packing:	100pcs/bag 2000pcs/ctn
Capability:	60,000pcs per day
Delivery time	7-30 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop