

Indication of [disposable](#) white PP non woven body-wrap:

Ideal for Travel, spa, beauty salon, beach, therapy center etc

The Specification of disposable white [PP non woven](#) body-wrap:

Name	Disposable white PP non woven body-wrap manufacturer
Material	PP, SMS
Style	Wrap on the chest, with or without elastic on the top
Size	S(80*100cm), M(80*110cm), L(100*125cm), XL(110*140cm),
Weight	20-65gsm is available(as your request)
Color	white, blue, green, pink, black, etc
Packaging	non sterile, 10pcs/bag, 50pcs/ctn or 100pcs/ctn,
MOQ	5000 PCS
Delivery	One 40HQ container can be finished within 15-25days
Store	Store in dry, humidity below 80%, ventilated, non-corrosive gases warehouse

The picture of disposable white PP non woven body-wrap:



orient Honest
MEDICAL INTERNATIONAL





Our service and advantage of disposable white PP non woven body-wrap:

1. Body-wrap Samples and customization available.
2. Professional QC to control the product quality. All the products comply with international quality standard such as CE,ISO13485, FDA.
3. Manufacturer with 15 years history, have enough experience and gained a high reputation from our customers all the world.
4. Quick payment, prompt delivery and competitive price.
5. Stable quality and best after-sale service.
6. OEM service is available.