

Disposable white non woven dustproof lab coat with pouch for factory

Lab coat applications and materials:

- 1) Hospital, Chemical and Pharmaceutical Industry
- 2) Electronic Room and Cleaning
- 3) Cleaning
- 4) Lab

Our **Lab Coat Regular Material:**

- 1) PP / SPP (100% Polypropylene Spun bond Non woven Fabric)
- 2) SMS (Polypropylene Spun bond Non-woven fabric + Melt blown Non-woven fabric + Polypropylene Spun bond (non-woven fabric))
- 3) PP + PE film
- 4) Micro

Lab coat design details

- 1) Collar: knit collar or turn down collar (single / double collar)
- 2) Front: buttons, plastic fastener cutout (plastic fastener piece), or front zipper available
- 3) Cuffs: open, elastic or knitted cuffs
- 4) Pocket: with / without pocket
- 5) Back situation: The back end can be opened about 20cm to make the user more comfortable

Lab coat size:

Coverall size	S	M	L	XL	2XL	3XL	4XL	5XL
Height (CM)	96	100	104	108	112	115	118	121
Chest circumference(CM)	140	144	148	152	156	160	164	168
Sleeve length(CM)	57	57	57	58	58	58	59	59

Lab coat packaging

Usually, 5 Pcs per bag, 50 Pcs or 100 Pcs per carton, or as your requirement.



[similar products:](#)

