

The property of disposable yellow **SMS isolation gown with thumb hole**

- The product is for single use only. After use please dispose according to local relevant laws and regulations
- The Flexible cuffs reduce chance of skin-debris being transferred to patient while maintaining dexterity
- It does not contain latex

The Specification of disposable SMS isolation **gown with thumb hole**

Fabric Material	SMS 45g
Style	open on back, cuff with thumb hole
Dimensions size	120*140cm, 130*150cm
Color	blue, yellow, or as your request
Certificate	CE, ISO13485, FDA
MOQ:	100,000pcs
Packing design	can be as your design
Packing	10pcs/bag, 1000pcs/ctn, or as your requirement
Capability	50,000pcs per day
Delivery time	5-15 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop

The picture of disposable SMS isolation gown AAMI Level 3 standard



