

Properties of EN14683 Type IIR **medical face mask** in blue

All masks are made of high quality material, providing excellent protection and superior comfort:

- Smooth and comfortable non-rolled material to eliminate skin irritation
- A very breathable construction layer stays cool, even during long cases
- Anti-fog masks are available with foam or film to prevent fogging
- Extra-long ties make donations easier

Specifications **EN14683 Type IIR** blue medical face mask

Dimension Size	17.5 * 9.5 cm
Fabric weight	PP25gsm + MB25gsm + PP25gsm, can be as your request
Color	Blue
Certificate	CE, ISO13485, FDA
MOQ:	100000 pcs
Packaging design	can be as your design
Packing	50 pcs / box, 2000 pcs / ctn, or as your requirements
Ability	1,000,000 pcs per day
Delivery time	7-30 days
Application	hospital, laboratory, clinic

Image of EN14683 Type IIR **medical face mask** in blue



