

Product details

Pattern of EO sterile sms surgical gown

- 1.Long sleeves/short sleeves/non sleeves
- 2.Elastic/knitted cuff
- 3.Ties on neck and waist
- 4.With sterilized and non-sterilized
- 5.Different specifications are available.
- 6.EO sterlie packing for hospital use.

Specification of EO sterile sms surgical gown

Material	SMS
Weight	50gsm
Style	Tie on / Velcro on the neck, elastic / knitted cuff, standard or reinforced
Size	120x130cm,130x150cm,150x160cm,150x190cm, 170x170cm
Color	White, green, blue, yellow, pink ,red, grey, black etc.
Packing	Non sterile,10 pcs/bag,10 bags/ctn. sterile,1 pcs/sterile bag,100 pcs/ctn
Packing design	All inner box and carton printing can be as your design
Carton size	58X28X48 cm
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures and so on

Product display



orient Honest
MEDICAL INTERNATIONAL






orient Honest
MEDICAL INTERNATIONAL



Our business process



The regular shipment way



FAQ

- 1, Free sample.
- 2, Small order acceptable;
- 3, L/C, T/T, D/P terms available;
4. FOB, CNF, CIF, EXW terms available;
- 5, Competitive and stable price because we are factory.
- 6, Qualified products. 3rd part product inspection is acceptable.
- 7, OEM service. More options on material, pattern, size and packing;
- 8, Timely delivery. The best lead time is 25-30 day for a 40HQ container.
- 9, Professional export service. Documentation, clearance, sea transportation, all in one.
- 10, Good after sale service. We take care complaints as we would like to be your long-term supplier!