

## Full color disposable non woven 4ply face mask

4-layers,inside and outside layers are coor PP non woven,middle layer is high standard 99 filter paper.With Ear-loop style.

Weight:25+25+25+25g (comfortable,without glass fibres,latex free.)

size:17.5x9.5cm(adult)

### advantages

1. Very low resistance to breathing,free from foreign smell,without any stimulating feeling.
2. Nose bar can be adaptable,comfortable,without glass fibres,latex free.
3. PFE(Filtration efficiency of non oily particles) $\geq$ 30%.
4. BFE(Bacteria filtration efficiency) $\geq$ 95%(99%).
5. Use ultrasonic welding technology, comfortable wearing.
6. white,blue,green,red,pink,yellow can be available

### packing

Non sterile,50pcs/box,2000pcs/ctn, carton size:52x38x30cm,

### indications

1. Sterile products are used for hospital clinical medical personnel and patients.
2. Non-sterile products can be used for dust-free workshop,hospital, laboratory, food industry, electronic manufactures and so on.
3. To prevent the spread of body fluids and blood splash content and isolate dust, particle, alcohol, blood, bacteria and virus invading.

Payment term:T/T,L/C,D/A,D/P/ESCROW

Delivery time:as you quantity(we can finish 1 millon pcs per day)

FOB port:Wuhan/Shanghai/Ningbo/Guangzhou

### disposable surgical face mask Specifications

Description	Full color disposable non woven 4ply face mask
Material	4ply,
Weight	25+25+25+25g (comfortable,without glass fibres,latex free.)
Style	Elastic ear-loop (latex free)
Size	17.5x9.5cm(adult)
Color	white,green,blue,pink,yellow,pink,orange,purple,red(all colors can be available)
packing	Non-sterile:50 pcs/box,2000pcs/ctn
packing design	All inner box and carton printing can be as your design
carton size	52x38x38cm

Applications	hospital clinical medical personnel and patients. dust-free workshop,, laboratory, food industry, electronic manufactures and so on
MOQ	5000 PCS
Sample	Can be supplied free for your quality checking within one week
Payment term	T/T,L/C,D/A,D/P/ES