

High quality [Disposable Non Woven Bed Sheets](#)

Disposable bed sheet styles: Flat, elastic end, square edge

Size: 80X190 + 20cm, 90X200 + 20cm, 160X190 + 20cm, 160X200 + 20cm, 180X200 + 20cm ...

Packing: 5 PCS / bag, 5 PCS / AG or 1 PCS / bag, 50 PCS / CTN.

Specification High quality [Disposable Non Woven Bed Sheets](#)

Fabric	PP, PP / PE, SMS, SMMS,
Style	Flat, elastic end, square tip, roll
Dimension size	80x190 + 20cm, 90x200 + 20cm, 160x190 + 20cm, 160x200 + 20cm, 180x200 + 20cm
Color	white, green, blue, yellow, pink etc.
Certificate	CE, ISO13485, FDA
MOQ:	5000 pcs
Packaging design	All inner bags and carton printing can be as per your design
Packing	10 pcs / box, 10 bags / ctn or 5 pcs / bag, 50 pcs / bag
Ability	Capacity: 10,000 pcs / day,
Delivery time	1 * 40HQ disposable bed sheet can be finished within 10-15 days;
Application	hospital, laboratory, clinic, food industry, dust free workshop
FOB port	Wuhan / Shanghai / Ningbo / Guangzhou
Payment method	T / T, PayPal, West Union, L / C, D / A, D / P, ESCROW,