

## Product details

### Specifications of medical face mask

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black

### Our funny face disposable surgical mask can achieve following quality standard

- 1) Conform to EN14683: 2005
- 2) Bacterial filtration efficiency (BFE%) > 99%.
- 3) TYPE IIR standard, resistance to penetration by synthetic blood-120mmHg.

### Packaging of pleated simple design mask

50pcs/box,40boxes/ctn,52x19x42cm;

1\*20 GP container can hold 1,400,000 Pcs of masks,700 Ctns;


1\*40 HC container can hold 3,450,000 Pcs of mask, 1725 Ctns.

Or as your requirements.

## Product display



## Our business process



## FAQ

Q1: .How do you solve quality issues if any?

A: Our factory will pay the highest attention to quality issues. We will deeply investigate the root cause of the problems and find the best solution to avoid any future problems. Meanwhile we will corresponding compensation for our customers.

Q2: What is your products expiration date?

A: Around 3-5 years.

Q3: Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE,ISO13485,FDA. REACH and NELSON report.