

The specifications of [SMS disposable pillow cover with edgefold](#)

Item name	SMS disposable pillow cover with edgefold
Material	SMS
Colour	White(as requested)
Gram weight	20g(as requested)
Type	disposable
Size	66X40cm(as requested)
Package	1PC/BAG,20PCS/MIDDLE BAG,200PCS/CTN(as requested)
Certificate	CE,ISO13485
Usage	Widely use in Beauty Salon,SPA,Beauty Center,Private Inn, Hospital
Feature	Disposable,breathable, comfortable, light weight and convenient
Application	Available in different weight, size and packing as requested;Customer's samples and specifications are always welcome

Other terms of SMS disposable pillow cover with edgefold

MOQ	10,000 pcs
Price	Negotiable, we will give you the best price
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

We can also produce:

1> non woven face mask, [dust mask](#), non woven bouffant cap, non woven clip cap,non woven doctor cap,

2> non woven surgical gown, non woven isolation gown, non woven lab, coverall coat,non woven apron, coverall, PE apron.

3> non woven sleeve cover, [PE sleeve cover](#)

4> [non woven shoe cover](#), PE shoe cover, CPE shoe cover, PE boot cover, non woven boot cover.

5> bed sheet, pillow case cover, non woven briefs

Details pictures of SMS disposable pillow cover with edgefold



orient Honest  
MEDICAL INTERNATIONAL





**orient Honest**  
MEDICAL INTERNATIONAL