

IIR surgical face masks , tie on medical mask , disposable 3ply face mask vendor

Performance of surgical face masks

- 1) To prevent the spread of body fluids and blood splash content.
- 2) Conform to EN14683 and FDA
- 3) Bacterial filtration efficiency (BFE%) > 99%.

Quick details of tie one medical mask

- 1) Layer: 3ply, PP non woven + high filtration paper + PP non woven
- 2) Size: 19.5X7.5cm
- 3) Color: white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black, we have 11 different colors totally.
- 4) Ear loop: tie on
- 5) Nose bar: Plastic or metal
- 6) Normal packing: 50 pcs/box, 2000 pcs/ctn, or as your requirement

The trade terms of surgical face masks

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	1,000,000 pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

Why choose us ?

- 1) Directly manufacturer for tie on medical mask, factory price
- 2) NELSON, ISO, CE certificates, BFE > 99%, type II / Type IIR standard.
- 3) OEM packing design, service in 24 hours.
- 4) Latex free, easy breath, eco friendly. Used in clinic, hospital, pharmacy, food processing, beauty salon, electronics industry, etc.

Picture of 3ply face mask

