

The Specification of COOL FEELING disposable face mask

Fabric Material	Layer 1: super soft spunlace Layer 2: high filtration paper Layer3: ES material
Style	Ear loop
Dimensions size	17.5*9.5 cm
Fabric weight	25gsm + 25gsm + 35gsm, could be as your request
Color	white, baby blue, baby pink, gray
Certificate	CE, ISO13485, FDA
MOQ:	100000pcs
Packing design	can be as your design
Packing	50pcs/box,2000pcs/ctn, or as your requirement
Capability	60,000pcs per day
Delivery time	7-30 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop

The picture of COOL FEELING disposable face mask




One more layer of protection,
one more peace of mind





orient Honest
MEDICAL INTERNATIONAL

0.32 cool sense of value



test item			
instant cool value	0.2-0.3	0.3-0.4	>0.4
somatosensory	not abvious	cool eithout ice	Have a cold feeling

* Data from the laboratory, different colors, fabric batches will have a certain impact on the value