

### **The property of level 2 disposable isolation gown ultrasonic seam**

These full length gowns are made from a fluid resistant, cool and comfortable material to provide both protection and comfort.

They are an inexpensive alternative to lab coats yet have key features such as waist ties and elastic cuffs.

The knee-length blue coats is unisex.

### **The Specification of level 2 disposable isolation gown ultrasonic seam**

Fabric Material	PP+PE, PP coated PE film
Style	With Velcro or ties on collar, ties on waist
Dimensions size	120*140cm, 130*150cm
Color	White, blue, yellow, or as your request
Certificate	CE, ISO13485, FDA
MOQ:	100,000pcs
Packing design	can be as your design
Packing	10pcs/bag, 1000pcs/ctn, or as your requirement
Capability	50,000pcs per day
Delivery time	5-15 days
Applications	hospital, laboratory, clinic, food industry, dust-free workshop

### **The picture of level 2 disposable isolation gown ultrasonic seam**



