

Product details

The feature of surgical sterile sheet

| | |
|-----------|--|
| Fabric | Polypropylene,SMS,PP+PE,or Microporous |
| Style | Fitted elastic design,Ideal for all sizes and styles of shelter cots |
| Function | Moisture proof |
| Size | 83" L x 35" W, or as your requirements |
| Packaging | Individually packaged |

About EO Sterilization for sterile hospital bed sheet , hygiene disposable sterile bed sheet

- 1) The sterilization indication label is in the each bag;
- 2) The bed sheet is sterilized by EO gas,and is packed in paper plastic bag;
- 3) The color of the sterilization indication label will change color if the bed sheet has been sterilized perfectly;
- 4) Normally, our packaging way is : 1Pcs per paper plastic bag, 50 Pcs per carton.

We can also produce as your requirements.

Product display



orient Honest
MEDICAL INTERNATIONAL

