

Our Advantage:

- 1.FDA, CE,ISO13485 certification.
- 2.One-stop service :excellent disposable medical consumable items.
- 3.Welcome any OEM requirements.
- 4.Qualified products,100% new brand material ,safe and sanitary.
- 5.Offered free samples .
- 6.Professional shipping service if necessary.
- 7.Full series after-sale service system.

Specifications:

Item	Medical PP Non Woven White Disposable Surgical Caps For Hospital
Material	PP 10g
Style	single elastic
Size	21inch
Color	white
Packing	Normally: 100pcs/bag, 20bags/ctn, 2000pcs/ctn. one 20ft container can load about 900 cartons(1,800,000pcs), one 40hq container can load about 2200cartons(4,400,000pcs).
Packing design	All inner box and carton printing can be as your design
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures etc.
MOQ	100,000PCS
Quality Control	1.We will not start to make the products until you confirm the sample. 2.The packing will be water-proof, damp-proof and sealed