

- **Suitable for people:** Doctor, Students, Patient, Cooker, Cleaner, office workers, etc
- **Suitable for the occasion:** Hospital, Food industry, Dust free workshop, School, Clean room, house etc
- **Product introduction:** the product is made of Non-woven fabric such as ultrasonic heating, tie or hook on neck, back waist with ties.
- **It is designed** under non woven material. Hot size in 115\*137cm which can fit most of users, can also support customize size and style to guarantees comfortable wearing and usage. Tight nice and economical bag pack is popular for consumers. Sterile Packing hot in medical.
- **This type of disposable surgical gown** has good air permeability, comfortable to wear, dust prevention Environment friendly. Custom produce can be adjustable to fit comfortable without glass fibre, hypoallergenic

#### Specifications

Description	China Hubei Wholesale CE ISO approved sterile one time use medical SMMS scrubs
Material	non woven, material gram weight from 18-60g Tell what you want, we will recommend one best to meet your require
Feature	Anti skid , Anti static, Hypoallergenic , Comfortable, High breathability, fiberglass free, latex free, friendly
Style	Tie wasit
Size	L115*W137cm, 120*140cm, 135*142cm etc (Length is size from Shoulder)
Color	Blue, white, yellow, green etc
Applications	Hospital clinical medical personnel and patients. dust-free workshop, laboratory, food industry, electronic manufactures and so on
Sample	Can be supplied free for your quality checking within one week
Payment term	T/T, L/C, D/A, D/P/ESCROW
Delivery time	1*40HQ can be finished around 35 days. ( 100,000 pcs per day)
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou







nt Honest  
CA INTERNATIONAL





Life at Honest  
INTERNATIONAL

