

[Disposable lab coat](#) is made from hydrophobic SMS, PP, PP+PE. CPE, Microporous material, Latex-free; Abrasion-resistant; Low lint;

Available color: Pink, white, green, blue, yellow, red, black etc. Material weight: 40-65GSM  
Lab coats size: S 96\*140\*57cm M 100\*144\*57 cm L 104\*148\*57 XL 108\*152\*58



Disposable lab coats products Advantages

1. With a high level of fluid repellency.
2. Latex-free; Abrasion-resistant; Completely impervious fluid. A good barrier to blood, body fluids and pathogens.
3. With standard and reinforced style. Reinforced style with extra protection reinforcement on sleeve & chest which can be completely impervious fluid and alcohol repellent.
4. Designed for optimal protection, comfortable wearing.
5. Capacity: 5000pcs/day, 1\*40HQ can be finished within 10-15 days;

Disposable lab coats indications

Disposable lab coats are widely used in medical safety protective gowns, aseptic workshop gowns, protective isolation,

Mining electronics factory, food factory farm animal husbandry biohazard and so on.

Disposable lab coats normal packing: 5 pcs/bag, 10 bags/ctn;



orient Honest  
MEDICAL INTERNATIONAL

