

## Product Details:

1. Material: Non woven polypropylene 100%
2. Size:80\*190+20cm,90\*200+20cm,160\*190+20cm,160\*200+20cm,180\*200+20cm or customized size
3. Elastic: 5mm along side width+30cm on length or customized size
4. Colour: Blue; White ; other size available
5. Packing: 10pcs/bag, 100pcs/carton

## Indications

- 1, Widely used for SPA, Hotel, Massage, Salon, Medical one time use
- 2, To prevent the spread of body fluids and blood splash content, and isolate the dust, particle, alcohol, blood, bacteria and virus invading.

## Specifications:

|          |   |
|----------|---|
| Material | Non woven polypropylene 100%  |
| Color    | White, blue   |
| Style    | elastics on width   |
| Size     | 80*190+20cm,90*200+20cm,160*190+20cm,160*200+20cm,180*200+20cm or customized size |
| packing  | 10 pcs/box, 100pcs/ctn  |
| Weight   | 25-65gsm or customized  |
| Package  | PE bag coated, then paper carton package or customized                            |
| Feature  | Durable, soft, disposable   |
| Usage    | SPA, Hotel, Massage, Salon, Medicinal, Hospital, Home                             |

## Our Services:

1. Provide the highest quality products

Professional QC team strictly control each step of the production until loading.

## 2. The fastest delivery time

If urgent, 1\*40HQ can be finished within 5-7days!

Normally, 1\*40HQ can be finished within 10-15 days!

Busy season, 1\*40HQ can be finished within 25-30 days!

## 3. The most reasonable price

Continuously invest in improving our skilled workers, products and processes  
to minimize the cost of each step!

## 4.The sweet and quick service

Our workers will be always so patiently making a response of your any questions within 5-10 minutes.

## 5.We can also help to source connected products.

Besides all kinds of disposable medical products, we can also help your source the connected products as your demand!

<https://www.china-disposable.com/products/Disposable-Fitted-Bed-Sheet-with-Elastic.html#.YLnHdPkzaUk>