

## **Nonwoven pleated face mask , modern surgical masks , nurse and patient protective mask**

### **The description of face mask**

- 1) Our face mask are hugely in demand for their softness, hygienic nature and reasonable rates.
- 2) Spun bond, SMS, Melt blown material used
- 3) Available 1ply, 2ply, 3 ply and 4ply
- 4) Varied color options : white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
- 5) Bulk and single sterile packaging offered.
- 6) Price as per quantity

### **Properties of modern surgical masks**

- 1) High filtration:3-ply pleated design offers,0.05 micron filtration at greater then 99%.
- 2) Slim aluminum wire fits comfortably over the nose,and maintain maximum filtration at the same time.
- 3) Tie and ear loop style: two pairs of polypropylene ties insure a comfortable fit; ear loop style attaches fast and with elastic loops.

### **The trade terms of modern surgical masks**

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW,Paypal
Supply ability	1,000,000 pcs/day,1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

### **Packaging of face mask**

50pcs/box,40boxes/ctn,52x19x42cm;  
1\*20 GP container can hold 1,400,000 Pcs of masks,700 Ctns;  
1\*40 HC container can hold 3,450,000 Pcs of mask, 1725 Ctns.  
Or as your requirements.

### **Picture of protective mask**



orient Honest  
MEDICAL INTERNATIONAL

