

PP beard cover masks , nonwoven disposable surgical beard mask , disposable beard net

Parameters of beard cover masks

- 1) Material: PP non woven
- 2) Color: White,blue,green....etc
- 3) Ear loop: elastic ear loop
- 4) Size: 19 Inch, 21 Inch

Features of disposable surgical beard mask

- 1) It's very low resistance to breathing,free from foreign smell,without any stimulating feeling ;
- 2) Elastic ear loop can be adjustable,comfortable wearing ;
- 3) Shelf Life: At least 3 years (upon receipt)

Packaging of beard net

Specification	Normal packing	Carton size
19 Inch	100 Pcs/PE bag ,1000 Pcs/carton	26*22*32 cm
21 Inch	100 Pcs/PE bag ,1000 Pcs/carton	52*23*26 cm

Our advantage

- 1) we are manufacturer for beard net with 12 years history, have large advantage in the price and technology;
- 2) Strong production power:1*40 HQ container can be finished in 15 days;
- 3) Promising sales team: service in 24 hours
- 4) High quality insurance: ISO9001, ISO13485, CE certificated, FDA Registered.

Picture of beard cover masks



orient Honest
MEDICAL INTERNATIONAL

