

Polypropylene face mask disposable , fashion designer surgical face masks , anti-dust filter mask

Indications of face mask disposable

- 1) Sterile face mask is for hospital, clinical, medical personnel and patients;
- 2) Non-sterile face mask can be used for dust-free workshop, food factory and so on.

Description of face mask disposable

- 1) 17.5*9.5cm, one size fits more;
- 2) Extra-soft ear loops eliminate pressure to the ears;
- 3) The inner layer made of soft facial tissue, no dye, no chemical, gentle to the skin, which is high filtration , BFE 99.5, PFE 99.7% at 0.1 micron;
- 4) Fluid resistant, latex free, fiber-glass free, odorless;
- 5) TYPE IIR quality standard, comfortable for medical professionals and healthcare personnel;
- 6) They are mainly used in hospitals, healthcare centers, dental clinics and pharmacies.

Packaging of designer surgical face masks

50pcs/box, 40 boxes/ctn, 52x19x42cm;
1*20 GP container can hold 1,400,000 Pcs of masks, 700 Ctns;
1*40 HC container can hold 3,450,000 Pcs of mask, 1725 Ctns.
Or as your requirements.

The trade terms of designer surgical face masks

MOQ	100,000 pcs
Price	Negotiable
Payment terms	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Supply ability	1,000,000 pcs/day, 1*40 HQ container can be finished in 15 days
Delivery time	Normally 30 days after confirmed the order
Sample	Sample is free, can be prepared well in 2 days

Picture of filter mask



orient Honest
MEDICAL INTERNATIONAL

