

Protective face mask 3 Ply with earloop ; anti-dust cleanroom face mask ; colored surgical masks

Description of face mask 3 Ply with earloop

- 1) 17.5*9.5cm, one size fits more;
- 2) Extra-soft ear loops eliminate pressure to the ears;
- 3) The inner layer made of soft facial tissue, no dye, no chemical, gentle to the skin, which is high filtration , BFE 99.5, PFE 99.7% at 0.1 micron;
- 4) Fluid resistant, latex free, fiber-glass free, odorless;
- 5) TYPE IIR quality standard, comfortable for medical professionals and healthcare personnel;
- 6) They are mainly used in hospitals, healthcare centers, dental clinics and pharmacies.

The specifications of cleanroom face mask

Material	PP non woven + high filtration paper + PP non woven
Fabric weight	Could be as your request
Ear loop	Tie on or ear loop
Size	17.5*9.5cm, or as your requirements
Color	white/blue/green/pink/yellow/orange/purple /lemon/Capri blue/burgundy/black
Packing	50 pcs/box,2000 pcs/ctn, or as your requirement
Packing design	all inner box and outer carton printing can be as your design
Applications	hospital,laboratory,food industry,dust-free workshop,electronic manufactures and so on

We produce following kinds of colored surgical masks

- 1) 3ply or 2 ply non woven face mask
- 2) 2 ply or 1ply paper face mask
- 3) 4 ply active carbon face mask

Picture of face mask 3 Ply with earloop



orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL

