

## Summary

Medical blankets that cannot be implanted are made of non -woven and polytel in it.

Provide warmth and comfort during the emergency of first aid

Widely used for the use of single patients

Perfect for emergency use, first aid equipment, emergency care centers, hospitals, and other medical care facilities

This material is latex free

## Material:

100% polyester inside, outside of PP Nonwoven

Features: disposable, waterproof

Type: Flat Sheet, Friendly-Co-Friendly Blanket

Form: rectangle

Pattern: striped

Use: hotel, disposable, hospital

is\_customized: Yes

Weight: 220-700G, 220-700G

Season: four seasons

Product Name: Disposable Medical Blanket

Keywords: disposable

Color: blue, white or adjusted

Fabric: 100% Polyets, Non -Woven PP

Use: airlines, hospitals, nursing homes

Logo: Receive adjusted logo

Size: Receive adjusted size size

Design: accept special designs

Package: 1 PC/POLI BAG

Fast detail

Place of origin: China

Material: polyester + non -woven

Features: folded, disposable, massage, waterproof

Technique: Non-Woven

Form: square

Use: warming up, patients move

is\_customized: Yes

Supply ability

Supply capability: 100000 pieces/pieces per week

Packaging & shipping

Port: FOB Wuhan/Shanghai/Ningbo

Waiting time:

Quantity (cut)	1 - 100000	> 100000
EST. Time (day)	20	To be negotiated

goods	score
Place of Origin	China
Ingredient	Polyester + non -woven
Feature	Folded, disposable, massage, waterproof
Technique	Non-Annoying
Form	Box
Use	Warming up patients and moves
is_customized	Yes



Used in: Nursing Home Emergency Hospital

can be packaged individually

Product description







