

## SMS blue disposable non woven patient transfer sheet suppliers

### Material features

SMS disposable stretcher sheets:

It is made from hydrophobic SMS/Spunlace material ,Latex-free; abrasion-resistant; low lint; with a high level of fluid repellency; a good barrier to blood, body fluids and pathogens. Used for clean, neat, convenient, safe.

Available color: blue.

Material weight:20-100gsm.

### Specifications:

Item	SMS blue disposable non woven patient transfer sheet suppliers
Material	SMS
Weight	20-110gms
Style	square end style (latex free)
Size	40"x60",40"x90",30"x72" ,60"x104" ,54"x88" ,31"x77", 33"x89" , 40"x48", 30"x84",40"x72",32"x72". Or as your demands.
Color	blue
Packing	5pcs/bag, 50pcs/bag or 10pcs/bag,100pcs/ctn
packing design	All inner box and carton printing can be as your design
carton size	Standard
Applications	Hospital for emergency bed and surgical bed
MOQ	500PCS
Sample	Can be supplied free for your quality checking within 2-4days
Payment term	T/T, L/C, D/A, D/P, Western Union, ESCROW, Paypal
Delivery time	1*40HQ can be finished within 10-15 days
FOB port	Wuhan/Shanghai/Ningbo/Guangzhou
Certificate	CE, ISO13485, FDA