

The property of Soft Touch Facial Towel

All cotton towel are made of the highest quality materials, providing excellent protection and superior comfort:

- 1.Size: 20 x 20 cm each piece or customized
2. Made of cotton material, soft and comfortable, clean and sanitary.
3. Convenient and lightweight , easy to carry out .
4. Can be wet and dry, are great for fragile or sensitive skin
5. Great water absorption, skin-friendly and comfortable.
6. Cotton tissue can instead of face towel and will be much cleaner to dry your face after washing.

The Specification of Soft Touch Facial Towel

Name	Soft Touch Facial Towel
Material	100% cotton;spunlace nonwoven
Style	Plain, mesh, pearl and customization patterns available
Size	15*20cm,20*20cm ,or other size can be customized
Color	white
Packing	1roll/portable bag ;100pcs/pack ,can be customized
Carton Box	All inner bag and carton printing can be as your design
Carton Size	52x38x30cm /52x19x30cm
Applications	personal care, makeup removing, house cleaning, travel use, camping, similar outdoor activities and business occasions.
MOQ	1000 bags
Sample	Can be supplied free for your quality checking within 3 working days
Payment way	T/T,PAYPAL,Westunion,L/C, D/A,D/P,ESCROW
Delivery time	as you quantity
FOB Port	Wuhan/Shanghai/Ningbo/Guangzhou

The picture of Soft Touch Facial Towel





orient Honest
MEDICAL INTERNATIONAL





orient Honest
MEDICAL INTERNATIONAL



Plant fibres




Safe enough



Elastic



Skin friendly




Fibrous porous
structure



The power
of touch

7.8inch







 Remove the bacteria



 Clean computer screen and accessories



 Clean cellphone screen



 Clean Automobile Interiors



orient Honest
MEDICAL INTERVENTIONS

Face Towel for Baby

Soft

