

Product details

The feature of non-woven bed sheet cover,hospital rubber bed sheets

| | |
|-----------|--|
| Fabric | Polypropylene,SMS,PP+PE,or Microporous |
| Style | Fitted elastic design,Ideal for all sizes and styles of shelter cots |
| Function | Moisture proof |
| Size | 83" L x 35" W, or as your requirements |
| Packaging | Individually packaged |

We also supply:

- 1) Isolation gowns,patient gowns,visit gowns,lab coats,coverall,scrub suits,kimono
- 2) doctor's caps,face mask,bed sheets,bed covers,pillow covers,sleeve covers,shoes covers, boot covers etc.

Why choose us ?

- 1) Direct manufacturer of bed sheet with elastic , factory price ;
- 2) 10,000 Pcs/day,we can finish 1*40 HQ container in about 15days;
- 3) CE ISO FDA approved , stable quality
- 4) Service in 24 hours, sample can be finished in 2 days, sample free.

Product display



orient Honest
MEDICAL INTERNATIONAL

The regular shipment way



FAQ

Q1: .How do you solve quality issues if any?

A: Our factory will pay the highest attention to quality issues. We will deeply investigate the root cause of the problems

and find the best solution to avoid any future problems.

Meanwhile we will corresponding compensation for our customers.

Q2: What is your products expiration date?

A: Around 3-5 years.

Q3: Do your company have any quality management system such as ISO, EN or CE certificate(s)?

A: Yes, we have CE,ISO13485,FDA. REACH and NELSON report.