

disposable EO Sterile Customized Surgery Cesarean section drape Pack set Used for Clinic

Product name	disposable EO Sterile Customized Surgery Cesarean section drape Pack set Used for Clinic
Application	clinic, hospital etc.
Material	non-woven, SMS, SM/PE, Bi-lam, Tri-lam, PE film, SS etc
Weight	30-50g sm etc
Product feature	1. Light weight and soft feeling
	2. Excellent water absorption and waterproof
	3. Blocking the transmission of bacteria
Size(cm)	See composition list or as request
Carton Size	60*40*45cm, or as request
Color	Blue
Certificate	CE & ISO13485
Sterilization	EO sterilized/non sterile
Payment term	T/T, L/C & D/P
Packing	1 set/blister bag, 10 pcs/carton, sterilized by EO, can be customized
Supply Ability	100000 pc/month
Delivery Time	20-40 days after order confirmed
Trade terms	FOB, C&F and EXW
Loading Port	Wuhan or Shanghai
OEM Project	Acceptable

Composition list				
No.	Item name	Description	Size	Unit
1	Drape utility with tape	45g SMS + adhesive tape	60*80cm	4
2	Medium drape	45g SMS	160*170cm	2
3	Back table cover	45g SMS + 55g PE film	115*225cm	1
4	Hand towel	66g white paper	40*40cm	2
5	Reinforced surgical gown	43g anti-static SMMS	120*150cm	2
6	C-section drape	45g SMS	195*275 cm	1
7	suture bag	45g transparent PE film	33*43cm	1
8	Mayo stand cover	PE film + 35g SMS	58*137cm	